

SOLUTIONS: Practice problems 2012

1. How would you prepare 400 ml of a 0.24 M NaCl solution (MW = 58.44 g/mole)?

$$[V \times C \times MW] \quad 0.4 \text{ L} \times 0.24 \text{ moles/L} \times 58.44 \text{ g/mole} = \underline{5.61 \text{ g NaCl dissolved in/brought up to 400 ml water}}$$

2. How would you prepare 750 ml of a 0.35 M Na₂PO₄ solution (MW = 141.96 g/mole)?

$$0.75 \text{ L} \times 0.35 \text{ moles/L} \times 141.96 \text{ g/mole} = \underline{37.265 \text{ g Na}_2\text{PO}_4 \text{ dissolved in/brought up to 750 ml water}}$$

3. How would you prepare 225 ml of a 0.15 M citric acid solution (453.6 g/mole)?

$$0.225 \text{ L} \times 0.15 \text{ moles/L} \times 453.6 \text{ g/mole} = \underline{15.309 \text{ g citric acid dissolved in/brought up to 225 ml water}}$$

4. What volume of stock 0.24 M NaCl is needed to prepare 75 ml of 0.1 M NaCl?

$$(c_1v_1=c_2v_2) \quad (0.24 \text{ M})(V_1) = (0.1 \text{ M})(75 \text{ ml}) \quad \rightarrow \quad V_1 = (0.1 \text{ M})(75 \text{ ml}) / (0.24 \text{ M}) = \underline{31.25 \text{ ml of stock NaCl}}$$

5. What volume of stock 0.35 M Na₂PO₄ is needed to prepare 150 ml of 0.2 M Na₂PO₄?

$$V_1 = (0.2 \text{ M})(150 \text{ ml}) / (0.35 \text{ M}) = \underline{85.714 \text{ ml of stock Na}_2\text{PO}_4}$$

6. What volume of stock 0.15 M citric acid is needed to prepare 230 ml of 2.5x10⁻³ M citric acid?

$$V_1 = (0.0025 \text{ M})(230 \text{ ml}) / (0.15 \text{ M}) = \underline{3.833 \text{ ml of stock citric acid}}$$

7. How much (ml) 3x10⁻³ M Na₂PO₄ can be prepared from 20 ml of 0.2 M Na₂PO₄?

$$[c_1v_1=c_2v_2] \text{ solve for } V_2: \quad V_2 = (0.2 \text{ M})(20 \text{ ml}) / 0.003 \text{ M} = \underline{1333.33 \text{ ml}}$$

8. If you have 35 ml of 0.1 M NaCl, is it enough to prepare 200 ml of 0.06 M NaCl?

$$[c_1v_1=c_2v_2] \text{ solve for } V_1: \quad V_1 = ((0.06 \text{ M})(200 \text{ ml}) / (0.1 \text{ M})) = \underline{120 \text{ ml}}, \text{ so no, you do not have enough stock.}, \text{ OR, using a different approach,}$$

$$\text{Solve for } V_2: \quad V_2 = ((0.1 \text{ M})(35 \text{ ml}) / (0.06 \text{ M})) = \underline{58.333 \text{ ml stock needed to make this amount.}}$$

9. For each of the following, indicate the mass (g) of reagent needed to prepare 125 ml of the indicated percent solution:

Reagent	Molecular weight	Percent solution needed (125 ml)	Mass of reagent (g)
NaCl	58.44	10 %	$(0.1 \text{ g/ml})(125 \text{ ml}) = 12.5 \text{ g}$
Na ₂ PO ₄	141.96	4.5 %	$(0.045 \text{ g/ml})(125 \text{ ml}) = 5.625 \text{ g}$
Citric acid	453.6	2.5 %	$(0.025 \text{ g/ml})(125 \text{ ml}) = 3.125 \text{ g}$

Percent concentration is a straight mass/volume (g/ml) calculation and doesn't need the MW. 1 g/ml = 100% solution.

10. You have 200 ml of 50X TAE. What volume of the 50X TAE is needed to prepare 25 ml of 10X TAE?

$$C_1V_1 = C_2V_2: \text{ solve for } V_1: V_1 = ((10X)(25 \text{ ml})/(50X)) = \underline{5 \text{ ml of 50X TAE stock}}$$

11. A protocol calls for using a working concentration of 1.5×10^{-4} M boric acid. To save shelf space, you decide to make up a 20X stock solution. What will be the molarity of that stock solution?

$$20 \times 1.5 \times 10^{-4} \text{ M} = \underline{3.0 \times 10^{-3} \text{ M}}$$

12. You make up 1000 ml of a buffer containing 300 ml citric acid and 700 ml disodium phosphate (Na₂PO₄). If the citric acid stock solution is 0.15 M and the Na₂PO₄ stock solution is 0.35 M, what are the working concentrations of these two components in the buffer?

$$\text{Citric acid} = (300 \text{ ml}/1000 \text{ ml}) \times (0.15 \text{ M}) = \underline{0.045 \text{ M or } 4.5 \times 10^{-2} \text{ M}}$$

$$\text{Na}_2\text{PO}_4 = (700 \text{ ml}/1000 \text{ ml}) \times (0.35 \text{ M}) = \underline{0.245 \text{ M or } 2.45 \times 10^{-1} \text{ M}}$$

13. You prepare a crude extract of alpha amylase from 2.7 g of germinated barley seed in 10 ml of buffer solution. What is the extract concentration in terms of milligrams of tissue extracted per milliliter of buffer?

$$2.7 \text{ g} / 10 \text{ ml} = \underline{0.27 \text{ g/ml or } 270 \text{ mg/ml}}$$

14. You use 0.5 ml of the extract prepared in # 13 in a total reaction volume of 6 ml. What is the working concentration of the extract in the reaction?

$$(0.5 \text{ ml})(270 \text{ mg/ml}) / 6 \text{ ml} = \underline{22.5 \text{ mg/ml}}$$

15. In a plate count analysis of *E. coli* cell density in a water sample you find 57 colonies on an LB agar plate inoculated with 100 ul of 10^{-4} diluted water sample. What is the cell density (#cells/ml) of *E. coli* in the water sample?

Given: 1 colony arises from 1 cell

$$(57 \text{ cells}/100 \text{ ul}) \times 10 = 570 \text{ cells/ml in the } 10^{-4} \text{ dilution} \quad (\text{this normalizes to cells/ml})$$

$$\underline{570 \text{ cells/ml} \times 10^4 = 5.7 \times 10^6 \text{ cells/ml} = \text{the original cell density in the water sample.}}$$